

Medik8®

SKINCARE BIOTECHNOLOGY™

Derma roller

New Titanium



Massively increase
product absorption.
Designed for home use.

Medik8® Titanium Dermaroller

Hailed as “the magic wand of skin care”, the new **Medik8 Titanium Dermaroller** is a fast and effective answer to all your skin concerns whether they are wrinkles, scars, acne scars, cellulite or stretch marks. Simple in notion and convenient for use, the **Medik8 Titanium Dermaroller** is designed to gently penetrate the top layers of your skin without causing damage or bleeding. It allows for much higher absorption of Medik8's topical formulation, leading to a significant increase in their efficacy. It also works with Medik8 serums to trigger the skin's natural healing response, leading to new collagen production and a firmer, tighter looking skin with diminished lines and wrinkles. **Medik8 Titanium Dermaroller** is constructed from pure titanium for a sharper edge and a virtually pain-free treatment experience. Take your roller home today.

Benefits

- No skin damage (epidermis is left intact)
- No skin sensitivity to UV light
- No risk of permanent skin discolouration
- Short healing time (no downtime)
- Painless and can be used at home
- Dramatically enhanced penetration of active ingredients
- Dramatically increased efficacy of topical skincare products


Results

When used in conjunction with Medik8 skincare

- Increased firmness and tightness of flaccid skin
- Thicker epidermis and enhanced dermal compactness
- Diminished wrinkles, lines and folds
- Improved skin texture and reduced pore size
- Reduced depth of indented acne or chicken pox scars
- Visibly relaxed old scars
- Lightened pigmentation marks and brighter skin tone
- Improved appearance of stretch marks and cellulite

Skin Needling Directions for Use:

1. Clean the Medik8 Titanium Dermaroller with Medik8 Sterilising Solution and allow to dry before every use.
2. Cleanse application area using Medik8 oneCleanse or your preferred Medik8 cleanser suitable for your skin type, and pat skin dry.
3. Visually divide the skin into sections (if on the face: forehead, left cheek, right cheek, nose, chin, left side of the neck etc.) Using a constant, gentle rolling motion roll up to 10 times back and forth over each skin section.



4. Apply your preferred Medik8 Prevent and/or Correct product, following the product instructions. Use the roller in the evening, as many topical formulations require time to absorb fully.
5. Finish with Hydr8 Night application.
6. Ensure the roller is cleaned thoroughly after each use, run under hot water and disinfect before allowing it to dry and returning to the case. Always allow the roller to cool before use.
7. Repeat several times per week, depending on skin tolerance. Ensure any irritation from prior treatments has fully subsided before repeat use. Gradually increase frequency of use, as your skin becomes adjusted to needling.

Do not use more than once a day.